

ZX1

1.360



Authorised for car client transport!
Tested according to ISO7176-19




MEYRA®



Characteristics of a fully welded wheelchair



Weight optimisation



Dynamic mobility traits



Modern design

Pure dynamism

COST EFFECTIVENESS

- for long-term care, with adjustment options to carry out different functions, which saves subsequent costs
- comes with a comprehensive range of standard equipment, which reduces the need for expensive customisations
- minimal wear and tear thanks to the integral construction of the frame and the quality of the components

THERAPEUTIC BENEFITS

- the optimal smoothness and lightweight design facilitate the independence of the user
- can be adjusted to suit various functions and individual needs of the user thanks to precise adjustments of back, seat and chassis
- comprehensive range of equipment for long-term users
- compact sitting, moving and transport thanks to optimised frame
- the compact version of the ZX1, which is reduced by 30 mm, fulfils the particular requirements for a compact sitting position (anterior frame angle of 85°)
- optimum mobility traits: sligher, sportier, easier to tilt and turn

APPLICATIONS

- for permanent and long-term care
- for people with paraplegia
- very well suited for various conditions, also for severely disabled people
- for users lacking the function of their lower extremities (lateral transfer, propelled with hand rims)
- for very active wheelchair users, handbikers and frequent travellers