

SMART^S

2.370




*Suitable for passenger transport.
Tested according to ISO 7176-19*




MEYRA®




Folding system optimised to reduce size



Operating requires very little force



Individually adjustable chassis



Unique side guard with height-adjustable armrests that can be operated with one hand

New Active Line

COST EFFECTIVENESS

- maximum part compatibility
- minimal wear and tear; repairs rarely necessary thanks to the integral construction of the frame and the quality of the components
- rapid, long-term supply of spare parts thanks to extensive network of specialist retailers

THERAPEUTIC VALUE

- stabilising and mobilising of users with neurological conditions possible thanks to optimal smoothness, good support when sitting and light weight
- chair can be easily wheeled right up to everyday objects thanks to the straight design of the front section and the swing-away, removable leg rests
- folding mechanism for independent loading into vehicle
- can be adjusted for various functions for the user, e.g. by retrofitting an electric drive
- comprehensive range of equipment for long-term users

FIELD OF USE

- for permanent and long-term care
- very well suited for various conditions, also for severely disabled people
- for users with the ability to stand and walk independently to some degree (transfer via the stand)
- for neurological conditions